

INFRARED FOREHEAD THERMOMETER

Non-contact Measurement "0" Risk of Cross-infection



LFR30B










SAFE



FAST



CONVENIENT

- 
SAFE Non-contact body temperature measurement avoids cross-infection
- 
FAST Imported probe from Germany, temperature displays in 1 second
- 
°C/°F Switch freely between celsius degree and fahrenheit degree, suitable for different user habits
- 
COLOR INDICATOR Three colored backlights indicate the different temperature ranges visually
 <math> < 37.5^{\circ}\text{C}</math> Green / - 
TWO MODES Two temperature measurement modes: Human body / Object
- 
DATA RECORD 99 sets of data can be stored for tracking the changes of body temperature
- 
POWER SAVING Auto shut down when there is no operation for 60 seconds

SPECIFICATIONS

Measuring Position	Forehead
Measuring Distance	0cm~5cm
Display	LCD Screen
Resolution	0.1°C / 0.1°F
Measuring Range	Human Body: 32.0°C ~ 43.0°C Object: 0°C ~ 100°C
Accuracy	For body: $\pm 0.3^{\circ}\text{C}$ (34°C~43°C) $\pm 0.4^{\circ}\text{C}$ (32°C~33.9°C)
Operating Temperature & Humidity	Temperature: 5°C~40°C Humidity: $\leq 85\%RH$
Storage Temperature & Humidity	Temperature: -20°C~55°C Humidity: $\leq 95\%RH$
Power Supply	2 X 1.5V AAA bateries
Service Life	5 years